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Literatur

- Knittel, Th., Lüthgens, K.,; Higher sensitivity for the detection of adenomas by a Hb/Hp bed site test using a novel sampling device (2013)
- Tao, S., Brenner; H.; Well adjusted qualitative immunochemical faecal occult blood tests could be a promising alternative for inexpensive, high-quality colorectal cancer screening; European Journal of Cancer Prevention 2013, 22:305-310
- Sieg A., Thoms C., Lüthgens K., John M.R., Schmid Gayk H.; Detection of colorectal neoplasm by the highly sensitive hemoglobin-haptoglobin complex in feces; Int J Colorectal Dis (1999) 14: 267-271
- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.
- and Mortality Weekly Report; (1995) 45 (5) : 107-110
- Hundt S, Haug U, Brenner H; American College of Physicians 2009; 150(3):162-169

FD Hb/Hp complex professionell rapid test

Intended Use

Our **FD Hb/Hp complex professionell** rapid test is based on a visual, immunochromatographic procedure for qualitative detection of human hemoglobin and hemoglobin-haptoglobin-complex in stool samples. The test is intended for professional in vitro diagnostic use in hospitals, medical laboratories and medical practices.

Medical Importance

The test supplies reliable evidence of occult blood in stool samples and thus supports the early detection of bowel cancer or adenomas. The simple and non-invasive test procedure allows an effective early screening (using colonoscopies for additional confirmation of diagnoses) and is essential for reducing bowel cancer rate.

Immunochromatography tests on hemoglobin and the Hb/Hp-complex in stool samples offers a variety of advantages in comparison to common tests (stool guaiac test). Most important is the fact that reactions based on human antigen antibody connection are highly specific and remain unaffected by nutritional components.

Tests based on Guaiac have a risk of false-positive results caused by animal hemoglobin/myoglobin or because of food containing pseudo-peroxidase. False-negative results can occur if the peroxidase is inhibited or blocked by vitamine C or other antioxidants.

Furthermore, human hemoglobin is detected in 100-fold lower concentration as in biochemical tests.

The **Hb/Hp-complex professionell** rapid test plays an important role in the recovery of hemoglobin from lysed erythrocytes and is less sensitive to digestion through proteolytic enzymes during the intestinal passage. That means the hemoglobin-haptoglobin-complex is detectable even after longer intestinal passages. Therefore it is also possible to detect adenomas and bowel cancer located in proximal intestinal areas. Additionally several studies showed the higher sensitivity for colorectal adenomas.

Our **FD Hb/Hp complex professionell** rapid test combines the detection of both free hemoglobin and Hb/Hp complex and provides sensitive results in colorectal cancer screening.

Principle of the Test

The immunological **FD Hb/Hp complex professionell** rapid test is based on an immunochromatographic method, in which both hemoglobin and hemoglobin/haptoglobin-complexes are specifically recognized through specific antibody reactions.

FD Hb/Hp complex professionell rapid test

For Hb:

The hemoglobin contained in the stool sample reacts with specific monoclonal antibodies that are attached to gold particles. This complex is distributed on the membrane and reaches the test band (T) on which the anti-hemoglobin is present.

For Hb/Hp:

The antibody attaches specifically to the hemoglobin/haptoglobin complex, where there is no cross-reactivity with hemoglobin.

In the case of a positive result, the corresponding molecules (hemoglobin or hemoglobin/haptoglobin complex) from the stool sample binds to the gold loaded antibodies on test band (T) whereby a pink-red coloration of the band becomes visible. If the result is negative, there are no molecules that can bind to the test band (T) and no coloration is shown. At this the coloration of control band (C) guarantees that sample application and migration is done correctly and that test is valid.

Batch specificity

Standard in-house validations make sure that it is possible to combine stool collecting tubes and test cassettes of different batches. This combination has no erroneous impact on test results. Batches are compatible and interchangeable with each other.

Contents

- 20 individually packed test cassettes
- 20 sets faecal sampling sets with sample-taking tubes, instructions on taking stool samples and stool collectors
- test instructions

Additional Required Materials

- Watch
- gloves

Warnings

- Valid for in vitro diagnostic and professional use only.
- Do not use the test after expiration date.
- Only open the test packaging if the test is actually going to be performed.
- All patient samples should be handled as contaminated.
- The extraction solution contains a small amount of sodium acid.
- Patients should follow the instructions for a sample collection exactly and should not take any sample during menstruation or in case of bleeding hemorrhoids. It is recommended to use enclosed stool container in order to avoid contact of stool with blood from urine samples or interfering components from the toilet water.


Storage and stability of reagents

The test can be stored at 2° to 30°C in a sealed bag until expiration date. The expiry date is mentioned on both the packaging of the rapid test cassette and on the label of the tube.

FD Hb/Hp complex professionell rapid test

FD Dripping stool sample tube

Contents:

- (1) Ejecting piston
- (2) Sampling rod
- (3) Intermediate section with wiper
- (4) Dripping Tube
- (5) Dripping tip
- (6) Closing cap

Dosing tip collects 40 mg stool sample

FD Hb/Hp complex professionell rapid test


Specimen Collection and Preparation

1. **Collect** the stool using the stool collector enclosed.


2. Use sample-taking tube and **detach sample rod** (2) (with corrugated cap) from tube by turning the cap.

3. **Stick** white sampling rod (2) into the **stool** and withdrawing it with a turning. Take special care that opening of sampling rod is completely filled with stool.

4. **Put** white sampling rod (2) completely in the intermediate section (3) and close it with a **quarter turn in clockwise direction**. Take care that the tube is completely sealed.


FD Hb/Hp complex professionell rapid test


Procedure of the Test


1. First **turn** the yellow ejection piston (1) **clockwise a quarter rotation** and then push it simultaneously down to transfer stool sample from the sampling rod (2) into the extraction buffer.

2. **Shake** stool-collecting tube (ideally on vortexer) until stool specimen has completely **dissolved** into extraction buffer.

3. **Remove** grey closing cap (6).

4. **Discard** the first drop of the sample and add **three drops** of solution via dripping tip (5) into the round sample window on the test cassette. First results can be obtained after **8 minutes**. Strongly positive samples may be readable earlier. **Viscous stool samples** may lead to extended running time (up to **15 minutes**).




FD Hb/Hp complex professionell rapid test


Interpretation of results



Positive: 2 pink-colored lines appear; one in the control region (C) and one in the test region (T). A similar intensity of both lines means that the test is positive. A weak staining of (T) indicates a weak positive result. (T) can be strongly colored than (C), by highly positives samples. If the Test is positive, a further examination with different clinical methods is recommended.

Negative: Just two pink-colored lines appear in the control region (C). In the test region (T) no line appears.

Invalid: If no line appears in control region (C), it is a sign that the test is not running properly or that the test materials are not in Order. Too many solid particles in the stool sample can also lead to an invalid result. In this case repeat the test with a new test cassette or contact the manufacturer for technical support.

FD Hb/Hp complex professionell rapid test

Restrictions

The test is used for qualitative determination of human hemoglobin and hemoglobin-haptoglobin complex in stool samples.

Stool samples should not be taken during or 3 days before or after menstruation, in the case of bleeding caused by constipation, bleeding hemorrhoids, or in the case of medicine being administered rectally. False positive results could be obtained.

Blood-thinning medications, such as for example, acetylsalicylic acid or coumarin, as well as iron preparations can lead to bleeding which is not caused by tumors.

No special diet is necessary before the performance of the test.

A positive test occurs if human hemoglobin or hemoglobin-haptoglobin complex is present in the sample. In addition to intestinal bleeding, blood in stool may also have other origins, such as hemorrhoids, menstruation bleeding or the mixing of the stool sample with urine containing blood. Not all intestinal bleedings are the result of benign or malignant adenoids. The results obtained by performing this test should be verified by using other methods.

Since carcinomas and adenoids can bleed to different extents and intermittently it is also advisable to check several stool samples of the same patient.

Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended to use the stool container included in the test.

Performance Characteristics

Analytical sensitivity:

The limit of detection of the test are 60ng for free and 25ng for complex hemoglobin per milliliter of buffer solution.

Prozone effect:

The **FD Hb/Hp complex professionell** rapid test also works reliably at extremely high hemoglobin values (>500,000ng/mL).

Specificity:

The **FD Hb/Hp complex professionell** rapid test is specific for human hemoglobin or the hemoglobin-haptoglobin complex and shows no cross-reactivity with hemoglobin from cattle, pigs, rabbits, horses, or sheep at a concentration of up to 0.5 mg/mL in the extraction buffer. The hemoglobin-haptoglobin complex also displays no cross-reactivity with bilirubin, vitamin C, and horseradish peroxidase.

FD Hb/Hp complex professionell rapid test

Literature

- Knittel, Th., Lüthgens, K.,; Higher sensitivity for the detection of adenomas by a Hb/Hp bed site test using a novel sampling device (2013)
- Tao, S., Brenner; H.; Well adjusted qualitative immunochemical faecal occult blood tests could be a promising alternative for inexpensive, high-quality colorectal cancer screening; European Journal of Cancer Prevention 2013, 22:305-310
- Sieg A., Thoms C., Lüthgens K., John M.R., Schmid Gayk H.; Detection of colorectal neoplasm by the highly sensitive hemoglobin-haptoglobin complex in feces; Int J Colorectal Dis (1999) 14: 267-271
- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.
- and Mortality Weekly Report; (1995) 45 (5) : 107-110
- Hundt S, Haug U, Brenner H; American College of Physicians 2009; 150(3):162-169)

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SYMBOL EXPLANATIONS

IVD	Only for in-vitro diagnostic use		Only for single use
Cont.	Content		Expiry date
LOT	Lot number		Storage temperature

FROST DIAGNOSTIKA GmbH

Speyerer Straße 74

D-67166 Otterstadt

Telefon: +49 (0) 6232 600487 0

Telefax: +49 (0) 6232 600487 60

www.frostdiagnostika.de

eMail: info@frostdiagnostika.de